

# Distortion Check

Automated Analysis of Distortion in MRgRT


Model 603S

*The fast, user friendly, cloud based solution.*



## distortion check

SOFTWARE FOR EVALUATION OF IMAGE DISTORTION



Distortion Check is a new, cloud-based solution designed to quickly and automatically quantify distortion in MRI images. Used in conjunction with CIRS MRI Grid phantoms, the software provides the capability to quickly and accurately measure distortion throughout the entire image volume.

After detecting all grid intersections, the software registers either a CAD or CT scan ground truth to these MR detected control points. An interpolation is then performed to generate the 3D distortion vector fields.

Results can be reported in a variety of output formats including scatter plots, contour plots, box and whisker plots for trending and DICOM overlays that can be imported to TPS or other 3rd party software. The software algorithms will work with any grid configuration and CIRS employs proprietary 3D printing techniques that enable easy modification of grid phantoms to meet customer requirements.

### Features

- Quickly & automatically analyze complete MR data sets
- Unique grid phantoms provide physical control points throughout entire 3D image volume
- Density of control points optimized to bring interpolation close to linearity
- CIRS materials simulate distortion due to susceptibility and chemical shifts typical to clinical patient scans
- Cloud based solution frees user of operating system and hardware constraints
- Online deployment facilitates collaboration and easy review and portability of results

## CIRS

900 Asbury Ave • Norfolk, Virginia 23513 • USA • Tel: 757-855-2765 • [WWW.CIRSINC.COM](http://WWW.CIRSINC.COM)

# Large Field MR Image Distortion Phantom


Assess MR image distortion in a tissue like medium

The Model 604-GS provides a tissue-like imaging environment for assessing geometric distortion in whole-body MRI scanners. The 604-GS is designed to assess image distortion caused by  $B_0$  inhomogeneity and nonlinearity of the magnetic gradients. Unlike other air-filled phantoms, the liquid-filled 604-GS is also sensitive to chemical shifts and susceptibility artifacts, additional causes of distortion found when encountering density differences in diagnostic MRI and radiation therapy treatment planning.

The phantom features a unique orthogonal 3D grid that provides complete geometric data throughout the imaging volume. Users can fill the phantom with the AAPM TG 100 recommended solution, or with specialized solutions for high-field MRI or unique scan sequences. Ground-truth measurements of grid locations can be obtained through CT scanning of the phantom.

Patient induced magnetic inhomogeneity occurs when a patient is placed in the scanner due to magnetic susceptibility, which causes the tendency of a material to magnetize in the presence of a magnetic field. Regions of abrupt change in tissue density or voids between tissue and air are prone to high magnetic susceptibility which disturbs the magnetic field. Distortions due to susceptibility and smaller distortions due to chemical shifts are better qualified by a phantom that presents such distortions caused throughout the entire field of view (FOV).


The phantom is marked for ease of alignment to position lasers and is designed for use with both curved and flat gantry tables.



**MODEL 604-GS**

## Features

- With a tissue-like MRI signal, assess image distortion due to:
  - $B_0$  field inhomogeneities
  - Gradient nonlinearity
  - Chemical shifts
  - Susceptibility artifacts
- Compatible with Distortion Check software (Free trial. See back for details)
- Submillimeter control point detection
- Optimized 3D grid spacing fills entire volume
- 2152 control points
- Easy integration and laser alignment
- Ground truth files for grid available



## MODEL 604-GS INCLUDES

| QTY | COMPONENT DESCRIPTION   |
|-----|---|
| 1   | Large Field MR Image Distortion Phantom   |
| 1   | Filling Funnel  |
| 10  | Extra plastic screws  |
| 1   | Harness   |
| 1   | Plastic tray  |
| 1   | Complimentary 90 day license for 5 successful analyzed scans using Distortion Check Software: For instructions on how to create your account, go to <a href="https://www.cirsinc.com/software/distortion-check/">https://www.cirsinc.com/software/distortion-check/</a> |
| 2   | User Guides (604-GS & 603S Distortion Check)  |
| 1   | Foam-Lined Carry Case   |
| -   | 60-Month Warranty   |

|                           |                                       |
|---------------------------|---------------------------------------|
| <b>OVERALL DIMENSIONS</b> | 300 mm (L), 276 mm (H) , 330 mm (DIA) |
| <b>WEIGHT (DRY)</b>       | 17 lb. (7.7 kg)                       |
| <b>WEIGHT (FILLED)</b>    | 62 lb. (28.1 kg)                      |



## SOFTWARE OPTIONS

| PART NO. | DESCRIPTION  |
|----------|--|
| 603S-25  | Distortion Check software (license for 2 years or 25 successfully analyzed scans, whichever expires first) Phantom Serial Number Required  |
| 603S-50  | Distortion Check software (license for 2 years or 50 successfully analyzed scans, whichever expires first) Phantom Serial Number Required  |
| 603S-100 | Distortion Check software (license for 2 years or 100 successfully analyzed scans, whichever expires first) Phantom Serial Number Required |
| 603S-200 | Distortion Check software (license for 2 years or 200 successfully analyzed scans, whichever expires first) Phantom Serial Number Required |
| 603S-500 | Distortion Check software (license for 2 years or 500 successfully analyzed scans, whichever expires first) Phantom Serial Number Required |

# MR Image Distortion & Image Fusion **Head** Phantom

## Assess MR distortion & Image Fusion in Stereotactic Radiosurgery Planning

CIRS Model 603-GS was designed to assess MR image distortion in Stereotactic Radiosurgery Planning. It's also a useful tool for verifying image fusion and deformable image registration algorithms used in various treatment planning systems. The tissue equivalent, anthropomorphic design closely matches a clinical imaging scenario. The phantom can be imaged using X-ray, Computed Tomography and Magnetic Resonance. It images well with all MRI sequences tested to date, including T1 weighted, T2 weighted, 3D Time of Flight, MPRAGE and CISS.





The entire inter-cranial portion of the skull volume is filled with an orthogonal 3D grid of 2.5mm diameter cross-like shaped rods spaced 10mm (I-S), 10.5mm (AP) and 11mm (L-R). Extra material added in the grid intersections increase grid signal. Five extended axis-rods intersect at the reference origin of the grid. The end of each extended axis is fitted with CT/MR markers allowing for accurate positioning with lasers and co-registration of CT and MR image sets.



**MODEL 603-GS**

### Features

- Works with Distortion Check Software
- Provides a realistic anthropomorphic scenario for CT and MR imaging
- Presents simulated bony anatomy as rigid landmarks for image fusion
- Unique inter-cranial 3D grid design used to assess spatial distortion
- Special pads compatible with all fixation frames
- CT/MR markers facilitate positioning and image registration
- 859 Control Points



### MODEL 603-GS INCLUDES

| MODEL  | QTY | COMPONENT DESCRIPTION   |
|--------|-----|---|
| -      | 1   | MR Distortion & Image Fusion Head Phantom   |
| -      | 1   | ABS Cradle  |
| -      | 1   | Complimentary 90 day license for 5 successful analyzed scans using Distortion Check Software: For instructions on how to create your account, go to <a href="https://www.cirsinc.com/software/distortion-check/">https://www.cirsinc.com/software/distortion-check/</a> |
| -      | 1   | Custom Carry Case   |
| -      | 1   | User Guide  |
| -      | -   | 60-Month Warranty   |
| 038-20 | 1   | SRS Frame Support Cups (set of 4)   |

|                           |                       |
|---------------------------|-----------------------|
| <b>OVERALL DIMENSIONS</b> | 32 cm x 24 cm x 18 cm |
| <b>WEIGHT</b>             | 12 lbs (5.5 kg)       |



### SOFTWARE OPTIONS

| PART NO. | DESCRIPTION  |
|----------|--|
| 603S-25  | Distortion Check software (license for 2 years or 25 successfully analyzed scans, whichever expires first) Phantom Serial Number Required  |
| 603S-50  | Distortion Check software (license for 2 years or 50 successfully analyzed scans, whichever expires first) Phantom Serial Number Required  |
| 603S-100 | Distortion Check software (license for 2 years or 100 successfully analyzed scans, whichever expires first) Phantom Serial Number Required |
| 603S-200 | Distortion Check software (license for 2 years or 200 successfully analyzed scans, whichever expires first) Phantom Serial Number Required |
| 603S-500 | Distortion Check software (license for 2 years or 500 successfully analyzed scans, whichever expires first) Phantom Serial Number Required |

# Distortion Check Software

Automatically analyze complete MR data sets



Distortion Check is cloud-based solution designed to quickly and automatically quantify distortion in MRI images. Used in conjunction with CIRS MRI Grid phantoms, the software provides the capability to quickly and accurately measure distortion throughout the entire image volume.

## Features

- Simple, user friendly web interface
- Detect physical control points throughout the 3D image volume
- Web based pdf report in summary or detailed format to NEMA MS 12 standard recommendations.
- Output raw data or DICOM overlay files for use with 3rd party software
- Establish multiple user accounts with different permissions
- Easily analyze and track multiple machines, imaging sequences and phantoms
- Establish distortion tolerance thresholds specific to different imaging sequences
- Re-compute any scans acquired for different tolerance thresholds

