Mammographic Accreditation Phantom

Model 015

USER GUIDE

CIRS

900 Asbury Ave • Norfolk, Virginia 23513 • USA • Tel: 757-855-2765 • WWW.CIRSINC.COM

BACKGROUND INFORMATION


The Mammographic Phantom was designed to test, the performance of a mammographic system by a quantitative evaluation of the system's ability to image small structures similar to those found clinically. Objects within the phantom simulate calcifications, fibrous calcifications in ducts, and tumor masses. The Phantom is designed to determine if your mammographic system can detect small structures that are important in the early detection of breast cancer. Test objects within the phantom range in size from those that should be visible on any system to objects that will be difficult to see even on the best mammographic systems.

DESCRIPTION OF THE PHANTOM

The 4.4 cm thick Mammographic Phantom is made up of: a 7 mm wax block insert containing 16 sets of test objects, a 3.4 cm (approx. 1-3/8") thick acrylic base, and a 3 mm (1/8") thick cover. All of this together approximates a 4.2 cm compressed breast of average glandular/adipose composition. Included in the wax insert are aluminum oxide (Al_2O_3) specks to simulate micro-calcifications. Six different size nylon fibers simulate fibrous structures and five different size lens shaped masses simulate tumors.

TEST OBJECTS

| Fibers | |
|--------|--------------------------|
| 1. | 1.56 mm nylon fiber |
| 2. | 1.12 mm nylon fiber |
| 3. | 0.89 mm nylon fiber |
| 4. | 0.75 mm nylon fiber |
| 5. | 0.54 mm nylon fiber |
| 6. | 0.40 mm nylon fiber |
| Specks | |
| 7. | 0.54 mm Al_2O_3 speck |
| 8. | 0.40 mm Al_2O_3 speck |
| 9. | 0.32 mm Al_2O_3 speck |
| 10. | 0.24 mm Al_2O_3 speck |
| 11. | 0.16 mm Al_2O_3 speck |
| Masses | |
| 12. | 2.00 mm (thickness) mass |
| 13. | 1.00 mm (thickness) mass |
| 14. | 0.75 mm (thickness) mass |
| 15. | 0.50 mm (thickness) mass |
| 16. | 0.25 mm (thickness) mass |

PRODUCT IDENTIFICATION SYMBOLS

Manufacturer


Consult Instructions for Use

USE OF THE PHANTOM

- A. For Accreditation Program: Follow the detailed instructions furnished separately by the American College of Radiology (ACR).
- B. For general system monitoring: Place the phantom on the image receptor surface in the same position as a breast. The nipple indent marker should be positioned away from the chest wall, just as the nipple of the patient's breast would be positioned. Position the x-ray tube and compression device as you would for a craniocaudal examination of a patient's breast. Choose the kVp and mAs factors as you would use for an average 4.5 cm breast and make an exposure. If the phototimer is used, set the phototimer to an inward position under the wax portion of the phantom. Use the phototimer settings as you would for an average 4.5 cm compressed breast. Process the exposed film or xeroradiographic plate. This image is over or under exposed, make a suitable adjustment in your factors and repeat the exposure. This is an indication that adjustments may be necessary for patient imaging of these compressed breast thicknesses and should be checked.
- C. Interpretation of image. Once an acceptable image has been obtained, have one or more viewers identify and record which objects are seen. Use a magnifying glass to look for the structures as you would on a patient's mammogram. A good imaging system should be able to see the 3rd speck group, 4th fiber, and the 3rd mass.

CARE OF THE PHANTOM

The phantom is fragile to impact and is adversely affected by extended exposure to temperatures above 110°F. The phantom should be kept clean and stored in a cool dry place when not in use. Clean with a damp cloth when needed.

| Reordering Information | |
|------------------------|--|
| Part No. | Description |
| 015-01 | Replacement Acrylic Housing/Contrast Disc for Model Accreditation Phantom |
| 015-02 | Replacement Acrylic Top for Model 015 Mammographic Accreditation Phantom |
| 015-03 | Replacement Set of Six (6) Acrylic Contrast Discs for Model 015 Mammographic Accreditation Phantom |
| 015-1 | Replacement ACR Wax Insert Only for Model 015 Mammographic Accreditation Phantom |

WARRANTY

All standard CIRS products and accessories are warranted by CIRS against defects in material and workmanship for a period as specified below. During the warranty period, the manufacturer will repair or, at its option, replace, at no charge, a product containing such defect provided it is returned, transportation prepaid, to the manufacturer. Products repaired in warranty will be returned transportation prepaid.

There are no warranties, expressed or implied, including without limitation any implied warranty of merchantability or fitness, which extend beyond the description on the face hereof. This expressed warranty excludes coverage of, and does not provide relief for, incidental or consequential damages of any kind or nature, including but not limited to loss of use, loss of sales or inconvenience. The exclusive remedy of the purchaser is limited to repair, recalibration, or replacement of the product at manufacturer's option.

This warranty does not apply if the product, as determined by the manufacturer, is defective because of normal wear, accident, misuse, or modification.

NON-WARRANTY SERVICE

If repairs or replacement not covered by this warranty are required, a repair estimate will be submitted for approval before proceeding with said repair or replacement.

RETURNS

If you are not satisfied with your purchase for any reason, please contact your local distributor prior to returning the product. Visit <https://www.cirsinc.com/distributors/> to find your local distributor. If you purchased your product direct through CIRS, call Customer Service at 800-617-1177, email rma@cirsinc.com, or fax an RMA request form to 757-857-0523. CIRS staff will attempt to remedy the issue via phone or email as soon as possible. If unable to correct the problem, a return material authorization (RMA) number will be issued. Non-standard or "customized" products may not be returned for refund or exchange unless such product is deemed by CIRS not to comply with documented order specifications. You must return the product to CIRS within 30 calendar days of the issuance of the RMA. All returns should be packed in the original cases and or packaging and must include any accessories, manuals and documentation that shipped with the product. The RMA number must be clearly indicated on the outside of each returned package. CIRS recommends that you use a carrier that offers shipment tracking for all returns and insure the full value of your package so that you are completely protected if the shipment is lost or damaged in transit. If you choose not to use a carrier that offers tracking or insure the product, you will be responsible for any loss or damage to the product during shipping. CIRS will not be responsible for lost or damaged return shipments. Return freight and insurance is to be pre-paid.

WITH RMA NUMBER, ITEMS MAY BE RETURNED TO:

CIRS
Receiving
900 Asbury Ave,
Norfolk, Virginia, 23513 USA

| PRODUCT | WARRANTY PERIOD |
|--|-----------------|
| Model 015 - Mammographic Accreditation Phantom | 60 Months |

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