

Lumbar Training Phantom

For interventional pain management

The CIRS Model 034 Lumbar Training Phantom provides a realistic puncture practice phantom for use with fluoroscopic image guidance.

Phantom contains realistic anatomy to facilitate eye/hand coordination in a training environment.

Phantom can also be imaged under CT, MR, and ultrasound.



Model 034



Features

- Anthropomorphic L-spine anatomy
- Self sealing puncture membrane
- Ligaments, spinal cord, discs, skin, and soft tissue have differing softness to permit trainee to "feel" the way to the injection site.

Practice

- Lumbar epidural
- Caudal epidural
- Facet and nerve blocks
- Sacroiliac joint injection
- Lumbar diskography

Model 034

ANATOMICAL DETAILS:

- Ilium
- Sacroiliac joint
- Sacrum plus L2 thru L5
- Superspinous ligament
- Interspinous ligament
- Ligamentum flavum
- Discs - Anulus fibrosis and nucleus pulposus

APPROXIMATE DIMENSIONS:

38 x 25 x 18 cm

APPROXIMATE WEIGHT:

19 lbs.

