

## ACR Digital Mammography Phantom Reference Image

This USB storage device contains a reference image of the pink wax insert for your phantom. It is to be used for location and orientation of all test objects within the phantom. It is not representative of the image quality one should expect when using the phantom for routine QA. The insert should be imaged in the acrylic phantom body, which, by design, produces additional attenuation and scatter. Image quality and test object detection may be degraded from what appears in this reference image.

To ensure quality, a CIRS trained inspector carefully reviews and approves each phantom image. In some cases phantoms may include debris or air bubbles in non-critical locations. These will not adversely impact use of the phantom. They will typically be annotated with an arrow on the image. Rest assured that your phantom meets specification. If however, you have any questions or concerns, please refer to the user guide for information on how to contact CIRS customer service.